

Aristocrat 888^{WNLS} SHOPPRIDER

Powerchair

Shoprider's Aristocrat 888^{WNLS} is all about freedom. Freedom to move where you want when you want. Manoeuvrability and affordability are the features offered by the Aristocrat line of power chairs. The full range of Shoprider power chairs meet the individual needs of the consumer. Perfect for both indoor and outdoor applications.

- > Shoprider 8 inch mid-wheel powerchair
- > Swing away joystick
- > On-board charger with retractable cord
- > 21AH batteries
- > Available colours: red or blue



888^{WNLS}

Specs	Weight capacity	Dimensions LxWxH	Total wt. w/bat.	Motor	Speed	Battery	Tire	24v chrgr	Driving dist.	Turning radius	Safe climb°
888 ^{WNLS}	113KG/250LBS	37"X 22.5"X 38"	64KG/141LBS	2 X .5HP	7KPH/4.4MPH	2 X 21AH	F6" M8" R5"	4AMP/HR	24KM/15M	432MM/17"	9°/15%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

www.eclipsemedical.com