

Aristocrat P424^L

SHOPPRIDER

Powerchair

Shoprider's Aristocrat P424^L is all about freedom. Freedom to move where you want when you want. Manoeuvrability and affordability are the features offered by the Aristocrat line of power chairs. The full range of Shoprider power chairs meet the individual needs of the consumer. Perfect for both indoor and outdoor applications.



- > Shopprider 14 inch mid-wheel powerchair
- > Full suspension
- > On-board 4A charger with 50AH batteries
- > Deluxe high back seat
- > Red, blue, black and champagne interchangeable panels included



P424^L

| Specs | Weight capacity | Dimensions LxWxH | Total wt. w/bat. | Motor | Speed | Battery | Tire | 24v chrgr | Driving dist. | Turning radius | Safe climb° |
|-------------------|-----------------|------------------|------------------|------------|-------------|----------|--------------|-----------|---------------|----------------|-------------|
| P424 ^L | 158KG/350LBS | 41"X 25"X 46" | 115KG/253LBS | 2 X 1.34HP | 9.6KPH/6MPH | 2 X 50AH | F6" M14" R5" | 4AMP/HR | 34KM/21M | 572MM/22.5" | 9°/15% |

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

www.eclipsemedical.com