

# Aristocrat P424<sup>M</sup>

**SHOPPRIDER**

## Powerchair

Shopperider's Aristocrat P424<sup>M</sup> is all about freedom. Freedom to move where you want when you want. Manoeuvrability and affordability are the features offered by the Aristocrat line of power chairs. The full range of Shopprider power chairs meet the individual needs of the consumer. Perfect for both indoor and outdoor applications.



- > Shopprider 10 inch mid-wheel powerchair
- > Full suspension
- > On-board 4A charger with 50AH batteries
- > Deluxe high back seat
- > Red, blue, black and champagne interchangeable panels included



### P424<sup>M</sup>

Specs	Weight capacity	Dimensions LxWxH	Total wt. w/bat.	Motor	Speed	Battery	Tire	24v chrgr	Driving dist.	Turning radius	Safe climb <sup>o</sup>
P424 <sup>M</sup>	136KG/300LBS	40"X 22.5"X 44"	100KG/220LBS	2 x .8HP	8KPH/5MPH	2 x 50AH	F6" M10" R5"	4AMP/HR	34KM/21M	490MM/18"	9°/15%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

[www.eclipsemedical.com](http://www.eclipsemedical.com)