

# Crown 888<sup>SL/SE</sup>

**SHOPPRIDER**

## Luxury scooter

Shoprider's Crown 888<sup>SL/SE</sup> provides feature rich luxury in personal mobility. Experience the outdoors in unprecedented style, comfort and performance. The Crown Series are the most popular scooters in the Shoprider Canada portfolio.

- > **Special Edition** Shoprider 4-wheel scooter
- > 1.2HP motor with 11.5 inch tires
- > Full suspension, anchor hooks
- > Delta tiller with full light package
- > 50AH batteries with 4A charger
- > Available colours: red, black or silver



### 888<sup>SL/SE</sup>

Specs	Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motor	Speed	Battery	Tire	24v chrgr	Driving dist.	Turning radius	Safe climb <sup>o</sup>
888 <sup>SL/SE</sup>	158KG/350LBS	48" x 23" x 43"	100KG/220LBS	1.2HP	12KPH/7.45MPH	2 x 50AH	4 x 11.5"	4AMP/HR	40KM/25M	1500MM/59"	9°/15%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

[www.eclipsemedical.com](http://www.eclipsemedical.com)