

## Cascade Pillow Back

The HAMILTON lift chair by Eclipse Medical is our classic pillow back model, offering an unparalleled level of support and comfort. The adjustable fill pillow back allows the individual to customize each chair to their unique liking. The Hamilton chair is a market leader in aesthetics, performance and quality.

- > Pillow top arms and fill adjustable pillow back makes the HAMILTON the obvious choice for comfort
- > Three position, full chaise, large side pockets with head and arm covers included
- > Heat and vibrating massage "option" available\*
- > Kiln dried hardwood and laminate joined in traditional wood-working techniques; lap joints, corner braces, dowel joints
- > Metal meets wood with counter bolstered inserts and machine bolts. Ensuring years of quiet operation
- > Premium lift mechanism by Liggett & Platt. Years of production quality with a strong lifetime warranty
- > Heavy duty OKIN lift actuators are the gold standard of motion furniture
- > Removable back (no tools required)
- > Standard integrated battery back-up

See Fabric page for fabric and colour options



### Specifications


<b>Sit</b>		<b>Stand</b>		<b>General</b>	
Overall Width	37"	Top of Back to Floor	60"	Style	Cascade Back
Overall Front to Back	34"	Seat Edge to Floor	22"	Position Range	3
Inside Arms	18.5"			Pockets (wide x tall)	27" x 10"
Seat Depth	20"	<b>Recline</b>		Headrest Cover	Yes
Back Height	28"	Back Angle	150°	Arm Covers	Yes
Seat Edge to Floor	20"	Chaise Angle	175°	Battery Backup	Yes
Back Edge to Floor	45"	Overall Front to Back	67"	Weight Capacity	375 lbs
Arm Height to Seat	7"	Edge of Footrest to Floor	22"	Warranty	8 YEARS (PRO-RATED)
		Distance off Wall	23"	Carton (L x W x H)	30" x 34" x 24.5"
				GW	135

\*Adjusting fill pillow not recommended for heat and massage model.

All Eclipse Lift Chairs are available in all fabric and colour options.

## Flock Fabric


Flock fabric offers a “velvety soft feel” with great durability and easy cleaning. All of this adds up to a high performance material that consumers continue to demand.



Deep Sea



Chocolate Sprinkle



Mountain Moss



Bordeaux

## Micro Suede

Since microfibres are very fine, they can be tightly woven into an extremely high quality fabric. The density of the fibres in the fabric leads to many desirable characteristics:

- soft and luxurious
- breathable and fade resistant
- easy to clean and durable




Vin Rouge



Cocoa Mousse



Rich Navy Blue



Desert Sage

## Faux Leather

The advantages of faux leather is that it remains soft and supple with virtually no maintenance. As well, it is moisture and stain resistant. In comparison to “real” leather, it does not crack or fade as easily.



Black



Brown