

Imperial 889^{XL}SBN

SHOPPRIDER

Premium scooter

Shopperider's Imperial 889^{XL}SBN is designed specifically for the person who wants to stand out from the crowd. These unique scooters offer performance and features second to none. Innovative yet practical to meet the needs of the consumer who wants it all.

- > Shoprider 4-wheel scooter
- > Delta tiller with 2-speed (hi-lo)
- > 2HP motor with 12.5 inch tires
- > Full suspension, full light package
- > 75AH batteries with 8A charger
- > Available colours: red or silver



889^{XL}SBN

Specs	Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motor	Speed	Battery	Tire	24v chgr	Driving dist.	Turning radius	Safe climb ^o
889 ^{XL} SBN	227KG/500LBS	56" x 28" x 50"	140KG/309LBS	2.0HP	13.5KPH/8.4MPH	2 x 75AH	4 x 12.5"	8AMP/HR	48KM/30M	1780MM/70"	12°/20%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

www.eclipsemedical.com