

Imperial 889^{XLSN}

SHOPPRIDER

Premium scooter

Shopperider's Imperial 889^{XLSN} is designed specifically for the person who wants to stand out from the crowd. These unique scooters offer performance and features second to none. Innovative yet practical to meet the needs of the consumer who wants it all.

- > Shoprider 4-wheel cabin scooter
- > Windshield wiper with fluid reservoir
- > 2HP motor with 12.5 inch tires
- > Full suspension with light package and delta tiller
- > 75AH batteries with 8A charger
- > Available colours: red or blue



889^{XLSN}

Specs	Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motor	Speed	Battery	Tire	24v chgr	Driving dist.	Turning radius	Safe climb°
889 ^{XLSN}	181KG/400LBS	63" x 29" x 64"	173KG/381LBS	2.0HP	13.5KPH/8.4MPH	2 x 75AH	4 x 12.5"	8AMP/HR	48KM/30M	1780MM/70"	12°/20%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

www.eclipsemedical.com