

Modern

The LEXINGTON lift chair by Eclipse Medical is our latest modern style chair, offering sleek lines and curves for the discerning individual. The LEXINGTON provides uncompromising quality, solid performance and home-enhancing opulent style all at a very affordable price point.

- > Three densities of foam ensures every part of the chair supports and cushions you properly
- > Three position, full chaise, dual side pockets with head and arm covers included
- > Heat and vibrating massage "option" available
- > Kiln dried hardwood and laminate joined in traditional wood-working techniques; lap joints, corner braces, dowel joints
- > Metal meets wood with counter bolstered inserts and machine bolts. Ensuring years of quiet operation
- > Premium lift mechanism by Liggett & Platt. Years of production quality with a strong lifetime warranty
- > Heavy duty OKIN lift actuators are the gold standard of motion furniture
- > Removable back (no tools required)
- > Standard integrated battery back-up

See Fabric page for fabric and colour options




Specifications

| Sit | | Stand | | General | |
|-----------------------|-----|---------------------------|------|-----------------------|---------------------|
| Overall Width | 34" | Top of Back to Floor | 61" | Style | Contemporary |
| Overall Front to Back | 34" | Seat Edge to Floor | 22" | Position Range | 3 |
| Inside Arms | 19" | | | Pockets (wide x tall) | 12" x 9" |
| Seat Depth | 20" | | | Headrest Cover | Yes |
| Back Height | 28" | Recline | | Arm Covers | Yes |
| Seat Edge to Floor | 20" | Back Angle | 150° | Battery Backup | Yes |
| Back Edge to Floor | 42" | Chaise Angle | 175° | Weight Capacity | 375 lbs |
| Arm Height to Seat | 7" | Overall Front to Back | 67" | Warranty | 8 YEARS (PRO-RATED) |
| | | Edge of Footrest to Floor | 20" | Carton (L x W x H) | 30" x 34" x 24.5" |
| | | Distance off Wall | 22" | GW | 132 |

All Eclipse Lift Chairs are available in all fabric and colour options.

Flock Fabric

Flock fabric offers a “velvety soft feel” with great durability and easy cleaning. All of this adds up to a high performance material that consumers continue to demand.



Deep Sea



Chocolate Sprinkle



Mountain Moss



Bordeaux

Micro Suede

Since microfibres are very fine, they can be tightly woven into an extremely high quality fabric. The density of the fibres in the fabric leads to many desirable characteristics:

- soft and luxurious
- breathable and fade resistant
- easy to clean and durable




Vin Rouge



Cocoa Mousse



Rich Navy Blue



Desert Sage

Faux Leather

The advantages of faux leather is that it remains soft and supple with virtually no maintenance. As well, it is moisture and stain resistant. In comparison to “real” leather, it does not crack or fade as easily.



Black



Brown