

Tri-Back

The MADISON lift chair by Eclipse Medical is our contemporary model, offering a luxuriously plush Tri-Back for exceptional lumbar support and comfort. The MADISON chair with its deluxe style, colours and fabrics will complement any room with a touch of elegance and warmth.

- > Pillow-top arms and fill adjustable head pillow brings this contemporary chair into the 21st Century
- > Three position, full chaise, large side pockets with head and arm covers included
- > Heat and vibrating massage "option" available
- > Kiln dried hardwood and laminate joined in traditional wood-working techniques; lap joints, corner braces, dowel joints
- > Metal meets wood with counter bolstered inserts and machine bolts. Ensuring years of quiet operation
- > Premium lift mechanism by Liggett & Platt. Years of production quality with a strong lifetime warranty
- > Heavy duty OKIN lift actuators are the gold standard of motion furniture
- > Removable back (no tools required)
- > Standard integrated battery back-up

See Fabric page for fabric and colour options




Specifications

Sit		Stand		General	
Overall Width	37"	Top of Back to Floor	59"	Style	Sculpted-lumbar
Overall Front to Back	34"	Seat Edge to Floor	22"	Position Range	3
Inside Arms	18.5"			Pockets (wide x tall)	27" x 10"
Seat Depth	20"			Headrest Cover	Yes
Back Height	27"	Recline		Arm Covers	Yes
Seat Edge to Floor	20"	Back Angle	150°	Battery Backup	Yes
Back Edge to Floor	44"	Chaise Angle	175°	Weight Capacity	375 lbs
Arm Height to Seat	7"	Overall Front to Back	66"	Warranty	8 YEARS (PRO-RATED)
		Edge of Footrest to Floor	21"	Carton (L x W x H)	30" x 34" x 24.5"
		Distance off Wall	22"	GW	132

All Eclipse Lift Chairs are available in all fabric and colour options.

Flock Fabric

Flock fabric offers a “velvety soft feel” with great durability and easy cleaning. All of this adds up to a high performance material that consumers continue to demand.



Deep Sea



Chocolate Sprinkle



Mountain Moss



Bordeaux

Micro Suede

Since microfibres are very fine, they can be tightly woven into an extremely high quality fabric. The density of the fibres in the fabric leads to many desirable characteristics:

- soft and luxurious
- breathable and fade resistant
- easy to clean and durable




Vin Rouge



Cocoa Mousse



Rich Navy Blue



Desert Sage

Faux Leather

The advantages of faux leather is that it remains soft and supple with virtually no maintenance. As well, it is moisture and stain resistant. In comparison to “real” leather, it does not crack or fade as easily.



Black



Brown