

Marquis 778^{EL}

SHOPPRIDER

Mid-size scooter

Shopperider's Marquis 778^{EL} is part of our classic line of scooters with the perfect blend of style and performance. All of these mid-size scooters provide the same great features with the option of three-wheel manoeuvrability (778^S model), four-wheel stability (888^{SLN} model) or a combination of both with the 778^{EL} model.

- > Shoprider 3/4-wheel scooter
- > 36AH batteries with 3A charger
- > Electric seat lift and full light package
- > C-type delta tiller
- > 4-wheel stability with 3-wheel manoeuvrability
- > Available colours: red, silver or blue



778^{EL}

Specs	Weight capacity	Dimensions LxWxH	Total wt. w/bat.	Motor	Speed	Battery	Tire	24v chrgr	Driving dist.	Turning radius	Safe climb ^o
778 ^{EL}	136KG/300LBS	43.5"x 22"x 41"	73KG/160LBS	0.84HP	8KPH/5MPH	2 x 36AH	F8" R10"	3AMP/HR	32KM/20M	1100MM/43"	8°/14%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

www.eclipsemedical.com