

Button Back

The ORION lift chair by Eclipse Medical is our traditional “neat and petite” 4-button pillow back. This chair is for individuals looking for an option where living space can be an issue. The ORION sets the standard for superb value without compromising on performance or aesthetics.

- > Neat and petite chair perfect for those tight living situations
- > 2.5 position, full chaise, large side pockets with head and arm covers
- > Heat and vibrating massage “option” available
- > Kiln dried hardwood and laminate joined in traditional wood-working techniques; lap joints, corner braces, dowel joints
- > Metal meets wood with counter bolstered inserts and machine bolts. Ensuring years of quiet operation
- > Premium lift mechanism by Liggett & Platt. Years of production quality with a strong lifetime warranty
- > Heavy duty OKIN lift actuators are the gold standard of motion furniture
- > Removable back (no tools required)
- > Standard integrated battery back-up



See Fabric page for fabric and colour options


Specifications

Sit		Stand		General	
Overall Width	31"	Top of Back to Floor	56"	Style	4-Button Back
Overall Front to Back	32"	Seat Edge to Floor	22"	Position Range	2.5
Inside Arms	19"			Pockets (wide x tall)	19" x 10"
Seat Depth	18"			Headrest Cover	Yes
Back Height	25"	Recline		Arm Covers	Yes
Seat Edge to Floor	19"	Back Angle	135°	Battery Backup	Yes
Back Edge to Floor	40"	Chaise Angle	175°	Weight Capacity	275 lbs
Arm Height to Seat	5"	Overall Front to Back	60"	Warranty	4 YEARS (PRO-RATED)
		Edge of Footrest to Floor	22"	Carton (L x W x H)	30" x 28" x 26"
		Distance off Wall	17"	GW	112

All Eclipse Lift Chairs are available in all fabric and colour options.

Flock Fabric

Flock fabric offers a “velvety soft feel” with great durability and easy cleaning. All of this adds up to a high performance material that consumers continue to demand.



Deep Sea



Chocolate Sprinkle



Mountain Moss



Bordeaux

Micro Suede

Since microfibres are very fine, they can be tightly woven into an extremely high quality fabric. The density of the fibres in the fabric leads to many desirable characteristics:

- soft and luxurious
- breathable and fade resistant
- easy to clean and durable




Vin Rouge



Cocoa Mousse



Rich Navy Blue



Desert Sage

Faux Leather

The advantages of faux leather is that it remains soft and supple with virtually no maintenance. As well, it is moisture and stain resistant. In comparison to “real” leather, it does not crack or fade as easily.



Black



Brown