

Regal GK-8³

SHOPPRIDER

Compact/portable scooter

Shopperider's Regal GK-8³ sets the standard in portable mobility with the world's first complete connector-less chassis design for ease of assembly and disassembly. It's as easy as 1-2-3.

- > Shoprider 3/4-wheel scooter
- > Lightweight and portable
- > Easy as 1-2-3 to disassemble
- > 12AH batteries with 2A charger
- > 8 inch flat free tires
- > Red, blue and champagne interchangeable panels included



GK-8³

Specs	Weight capacity	Dimensions LxWxH	Total wt. w/bat.	Motor	Speed	Battery	Tire	24v chrgr	Driving dist.	Turning radius	Safe climb ^o
GK8 ³	114KG/250LBS	38"X 21.5"X 35"	41KG/89LBS	0.5HP	6KPH/3.8MPH	2 x 12AH	F7" R8"	2AMP/HR	14KM/9M	838MM/33"	6°/10%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

www.eclipsemedical.com