

# Regal GK-9

**SHOPPRIDER**

## Compact/portable scooter

Shopperider's Regal GK-9 sets the standard in portable mobility with the world's first complete connector-less chassis design for ease of assembly and disassembly. It's as easy as 1-2-3.

- > Shoprider 4-wheel scooter
- > Lightweight and portable
- > Easy as 1-2-3 to disassemble
- > 21AH dual split battery pack with 3A charger
- > 9 inch flat free tires
- > Available colours: red



### GK-9

Specs	Weight capacity	Dimensions LxWxH	Total wt. w/bat.	Motor	Speed	Battery	Tire	24v chgr	Driving dist.	Turning radius	Safe climb°
GK9	136KG/300LBS	41.5"X 22"X 35"	56KG/123LBS	0.75HP	7KPH/4.4MPH	2 X 21AH	F9" R9"	3AMP/HR	24KM/15M	1150MM/45"	6°/10%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

[www.eclipsemedical.com](http://www.eclipsemedical.com)