

Bigfoot S846

TRAILMASTER

Full-size scooter

The Eclipse Trailmaster series allows individuals to experience the outdoors in feature rich design, style, performance & comfort. They also come standard with an LCD digital dashboard. This series is for those individuals who want to travel a little further, in style.

The full-size series are the most popular scooters in our portfolio.

- > Trailmaster Bigfoot 4-wheel scooter
- > Digital LCD dashboard
- > Full light package
- > Front & rear suspension w/ 400 lbs. weight capacity
- > 75AH batteries with 6A charger
- > Colours: red, blue, orange or green



\$5,995.00
MSRP



S846 Specifications

Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	Turning radius	Safe climb°
182KG/400LBS	56"x 26"x 49"	130KG/286LBS	2.3HP	15KPH/9.5MPH	2 x 75AH	F12" R12"	50KM/32M	161CM/63"	10°/16%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

www.eclipsemedical.com

#052025EN