

# Dart P321B

# PATHMASTER

## Powerchair

Eclipse Pathmaster Series – Freedom to move where you want when you want. Manoeuverability and affordability are the features offered by our line of power chairs. The full range of power chairs meet the individual needs of the consumer. Perfect for both indoor and outdoor applications.



- > Pathmaster Dart lightweight powerchair
- > Compact stylish design
- > 22AH batteries
- > Off-board 4A charger
- > Colours: red, blue or white



**\$3,395.00**  
MSRP

### P321B Specifications

Weight capacity	Dimensions LxWxH	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	Turning radius	Safe climb <sup>o</sup>
136KG/300LBS	37.5"x 22"x 36"	51.4KG/113.5LBS	2 x .4HP	6.4KPH/4MPH	2 x 22AH	F7" R8"	24KM/15M	40CM/15.5"	6°/10%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

[www.eclipsemedical.com](http://www.eclipsemedical.com)

#012023EN