

Interceptor S940GTX

INTERCEPTOR

Premium scooter

The Eclipse Interceptor series of premium style scooters takes performance, features and comfort to new levels. Modern chassis design and lighting packages allow the most discerning individuals to travel in ultimate style. Great for the outdoors to visit people and places near and far.

- > **New Digital LCD Dashboard Display** featuring:
Speedometer, Odometer, Temperature, Time and Battery Gauges
- > **Modern L.E.D. Lighting Package** includes:
Dual Head, Tail, Signal, Hazard and Brake lighting system
- > Robust drivetrain powered by a P&G 140 amp S-Drive controller with a high RPM four-pole motor. Delivers exceptional performance.
- > Delta tiller (wrap around handle bar). Fits two rear view mirrors. Comes equipped with an extra large glove box.
- > Colours: red, blue or silver



\$6,995.00
MSRP



* Batteries not included.

S940GTX Specifications

Weight capacity	Dimensions LxWxH	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	Turning radius	Safe climb ^o
204KG/450LBS	61"x 29"x 47"	165KG/363LBS	3HP	15KPH/9.5MPH	FITS 75-100	F14" R16"	40KM/25M	180CM/71"	10°/16%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

www.eclipsemedical.com

#012023EN