

Pioneer S148

TRAILMASTER

Mid-size scooter

The Eclipse Trailmaster series offer modern design, style, performance & features such as an LCD digital dashboard. This mid-size line of products are for individuals looking to access their local community for shopping, appointments and visiting with family and friends.

- > Trailmaster Pioneer 4-wheel scooter
- > Digital LCD dashboard
- > Light package with active brake lights
- > Rear suspension
- > 50AH batteries with 5A charger
- > Colours: red, blue or silver



\$4,995.00
MSRP



S148 Specifications

Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	Turning radius	Safe climb ^o
240KG/450LBS	51"x 25"x 45"	110KG/242LBS	2.3HP	12.8KPH/8MPH	2 x 50AH	12"	40KM/25M	148CM/57"	10°/16%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

www.eclipsemedical.com

#012023EN