

Spyder-R P326AR Rehab

PATHMASTER

Powerchair

Eclipse Pathmaster Series – Freedom to move where you want when you want. Manoeuverability and affordability are the features offered by our line of power chairs. The full range of power chairs meet the individual needs of the consumer. Perfect for both indoor and outdoor applications.

- > Pathmaster Spyder-R 10 inch mid-wheel powerchair
- > Full suspension
- > Off-board 4A charger with 36AH batteries
- > 16" – 17" – 18" width and depth adjustable
- > Tension adjustable back
- > Swing-away joystick
- > Colours: red or blue



\$5,990.00
MSRP



P326AR Rehab Specifications

Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	Turning radius	Safe climb ^o
136KG/300LBS	34"x 24"x 40"	88KG/194LBS	2 x .7HP	8KPH/5MPH	2 x 36AH	F6" M10" R6"	28KM/18M	50CM/20"	6°/10%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

www.eclipsemedical.com

#012023EN