
Axis UL-8WSLA



Powerchair

Freedom to move where you want when you want.

Manoeuverability and affordability are the features offered by our line of power chairs. The full range of power chairs meet



- Shoprider ultra-light powerchair
- Easy as 1-2-3 to disassemble
- 12AH batteries
- Off-board 2A charger
- Colour: red







UL-8WSLA Specifications

Total wt. w/bat. Safe climbo Weight capacity **Dimensions** LXWXH Motor **Speed Battery** Tire Driving dist. Turning radius 113kg/250LBS 30"x 21.5"x 32" 47kg/104LBS 2 x .5HP 6крн/3.9мрн 2 x 12ан F6" R9" 14км/9м 40cm/15.5" 8°/14%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again-in style