

Chameleon GK-8⁴

SHOPPRIDER

Compact/portable scooter

Setting the standard in portable mobility with the world's first complete connector-less chassis design for ease of assembly & disassembly. It's as easy as 1-2-3. This class of scooters are ideal for individuals going shorter distances and needing to transport their unit to get from place to place.

- > Shoprider 4-wheel scooter
- > Lightweight and portable
- > Easy as 1-2-3 to disassemble
- > 12AH batteries with 2A charger
- > 8 inch flat free tires
- > Colours: red, blue or champagne with interchangeable panels included



\$2,295.00
MSRP



GK-8⁴ Specifications

Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	Turning radius	Safe climb ^o
114KG/250LBS	39" x 21.5" x 35"	43KG/94LBS	0.5HP	6KPH/3.8MPH	2 x 12AH	F7" R8"	14KM/9M	94CM/37"	6°/10%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

www.eclipsemedical.com

#012023EN