

# Cobra 778<sup>EL</sup>

**SHOPPRIDER**

## Mid-size scooter

Our classic line of scooters with the perfect blend of style and performance – for individuals looking to access their local community for shopping, appointments and visiting with family and friends. All of these mid-size scooters provide the same great features with the option of 3-wheel manoeuvrability (778S), 4-wheel stability (888SLN) or a combination of both with the Shoprider 778EL model.

- > Shoprider 3/4-wheel scooter
- > C-type delta tiller
- > Electric seat lift and full light package
- > 36AH batteries with 5A charger
- > 4-wheel stability with 3-wheel manoeuvrability
- > Colours: red, silver or monarch blue



**\$3,795.00**  
MSRP



### 778<sup>EL</sup> Specifications

Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	Turning radius	Safe climb <sup>o</sup>
136KG/300LBS	43.5" x 22" x 41"	73KG/160LBS	0.84HP	8KPH/5MPH	2 x 36AH	F8" R10"	32KM/20M	110CM/43"	8°/14%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

[www.eclipsemedical.com](http://www.eclipsemedical.com)

#092019EN