

Navigator^L P424^L

SHOPPRIDER

Powerchair

Freedom to move where you want when you want.

Manoeuverability and affordability are the features offered by our line of power chairs. The full range of power chairs meet the individual needs of the consumer. Perfect for both indoor and outdoor applications.

- > Shoprider 14 inch mid-wheel powerchair
- > Full suspension
- > Off-board 5A charger with 50AH batteries
- > Swing-away joystick
- > Colours: red, blue, black or champagne with interchangeable panels included



\$4,995.00
MSRP



P424^L Specifications

Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	Turning radius	Safe climb ^o
158KG/350LBS	41" x 25" x 46"	115KG/253LBS	2 x 1.34HP	9.6KPH/6MPH	2 x 50AH	F6" M14" R5"	34KM/21M	57CM/22.5"	9°/15%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

www.eclipsemedical.com

#092019EN