

Pirouette 888^{WNLS}

SHOPPRIDER

Powerchair

Freedom to move where you want when you want.

Manoeuverability and affordability are the features offered by our line of power chairs. The full range of power chairs meet the individual needs of the consumer. Perfect for both indoor and outdoor applications.

- > Shoprider 8 inch mid-wheel powerchair
- > Swing away joystick
- > Off-board 5A charger with 21AH batteries
- > Colours: red or blue



\$4,295.00
MSRP



888^{WNLS} Specifications

Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	Turning radius	Safe climb ^o
113KG/250LBS	37" x 22.5" x 38"	64KG/141LBS	2 x .5HP	7KPH/4.4MPH	2 x 21AH	F6" M8" R5"	24KM/15M	43CM/17"	9°/15%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

www.eclipsemedical.com

#012023EN